

Hospital Bed

ABS



- Sides and floors are made from ABS
- Sides are Foldable
- Tilting/Inclination 0-90°
- Knee rest 0-45°
- Removable head / foot Panel (ABS)
- Bear Weights up to 200 Kg
- Wheels (Diameter:14 cm) with brakes
- Dimensions: (93*198.5*53.5)(width*Length*height)
- IV Pole

Office: 17 Ali Youssef st. - Elqasr Alainy, Cairo

Factory: New Borg Alarab - 2nd Industrial zone - Block 18

E-mail: info@agcfactory.com

sales@agcfactory.com

Mobile: 01274990092

Website: www.agcfactory.com