

# Phototherapy

Model: PL



Phototherapy is used for jaundice treatment in newborn infants who have hyperbilirubinemia

It can be used on newborn infant bed or above the incubator, which phototherapy head can be directed to infant body

- 5 LED Phototherapy Lamp (Blue Light)
- Wavelength for LED lamps is 460 nm
- Counter to count LED lamp life time and can be reset after each session
- Height can be adjusted
- Wheels can be locked
- LED lamp life time is 30000 hours
- Angel for Lamp module can be adjusted



Office: 17 Ali Youssef st. - Elqasr Alainy, Cairo

Factory: New Borg Alarab - 2nd Industrial zone - Block 18

E-mail: [info@agcfactory.com](mailto:info@agcfactory.com)

[sales@agcfactory.com](mailto:sales@agcfactory.com)

Mobile: 01274990092

Website: [www.agcfactory.com](http://www.agcfactory.com)