


# x-ray Machine



- **Accurate**

- **Affordable**

- **Reliable**


- **Simple Digital**

## Generator & Module



- **Best fit and cost effective for medical and veterinary applications**
- **Auto calibration for various tubes and cable length**
- **Suitable for Film, CR and DR applications**
- **Auto exposure control(AEC) function to calculate cumulative exposure and limit excessive radiation dose**
- **Compatible with over variable X-ray tube models using low speed starter**

# Generator & Module




- Auto calibration for various tubes and cable length
- APR settings over numerous techniques
- Automatic settings of optimal technical conditions for the best images
- Art anatomical graphics selection with touch screen

## Generator & Module

Generator type	Radiography
Max. Output	40kW
kV Range(Rad.)	40-125(150)kV
kV Range(Fluoro.)	-
mA Range(Rad.)	10-500mA
mA range(Fluoro.)	-
mAs Range	0.1-630mAs
Number of APR	216 (1096)
Power Supply	1Ø 220/240VAC
Exposure Time(Rad.)	0.004-10 sec
Exposure Time(Fluoro.)	-
Option	AEC
Peripheral Interface	S, P, E*


# Tube Head Features

## Balanced Motion



Clear Tube Head Panel  
Easy Tube Moving

# Wall Bucky Features




**Compatible with various selections of DR systems**

**Smooth Sliding Motions  
Up & Down**

- \*Contains a cassette cabinet that moves up and down from ground level up to 200 cm in sliding motion system**
- \*And moving weights, with mounting brakes (manual).**
- \*Variable cassette holder sizes from 8 x 10 inches up to 17 x 17 inches**
- \*Contains linear ray filter orthogonal lines 103 line ratio 1:10 made in Korea**

# 4way Table



- \* Imaging table with mechanical compass
- \* Moving table in 4 directions and weights up to 180 kg.
- \* Sizes of the cassette place from 8 x 10 inches up to 17 x 17 inches
- \* Fine grid x-ray filter orthogonal lines 103 lines ratio 1:10
- \* Installation of the imaging table does not require any installation ceiling or walls or any modification in the room

Office: 17 Ali Youssef st. - Elqasr Alainy, Cairo

Factory: New Borg Alarab - 2nd Industrial zone - Block 18

E-mail: [info@agcfactory.com](mailto:info@agcfactory.com)      [sales@agcfactory.com](mailto:sales@agcfactory.com)

Mobile: 01274990092 - 01002219714

Website: [www.agcfactory.com](http://www.agcfactory.com)

