

AGC



AGC
For
Technological
Industries



MR-200
X-Ray Machine
200 mA

Office : 17 Ali Youssef st. Elqasr Alainy-Cairo, Egypt

Factory : Alex-New Borg Alarab-2nd Industrial Zone-Block 18

info@agcfactory.com

sales@agcfactory.com

+201007700679

+201141199418

www.agcfactory.com



AGC

AGC For Technological Industries

Control

- * 100 KVP - 200 mA
- * Normal frequency 220 V — 50 HZ
- * Contains the keys to control the following factors :
 - * KV : From 40 to 100 kVP. Control up to 1 kVP
 - * mA : from 50 to 200 mA
 - * Time : From 0.01 to 3 seconds

X-Ray Tube

- * Jx51-20, 40/125 made in China 200 mA

Tube Stand

- * The tube holder moves vertically to change the shooting distance and moves horizontally on the longitudinal axis of the table

4 Way Floating Table

- * Imaging table with mechanical compass
- * Moving table in 4 directions and weights up to 180 kg.
- * Sizes of the cassette place from 8 x 10 inches up to 17 x 17 inches
- * Fine grid x-ray filter orthogonal lines 103 lines ratio 1:10
- * Installation of the imaging table does not require any installation ceiling or walls or any modification in the room

Chest Stand

- * Contains a cassette cabinet that moves up and down from ground level up to 200 cm in sliding motion system
- * And moving weights, with mounting brakes (manual).
- * Variable cassette holder sizes from 8 x 10 inches up to 17 x 17 inches
- * Contains linear ray filter orthogonal lines 103 line ratio 1:10 made in Korea